

| Detector Position | Laser | Filter | Example Fluorophores |
|-------------------|--------|--------|---|
| FL1 | 488 nm | 530/30 | FITC, Alexa Fluor® 488, CFSE, GFP, YFP |
| FL2 | 488 nm | 585/40 | PE, PI, PE-CF594 |
| FL3 | 488 nm | 670 LP | PerCP, PE-Cy5, PerCP-Cy5.5, PE-Cy7, 7-AAD |
| FL4 | 640 nm | 675/25 | APC, Alexa Fluor® 647 |